

**MRT-PCR Purification using Exonuclease***Influenza Division, Virology, Surveillance and Diagnosis Branch, Genomics and Diagnostic Team*

**NOTE: This procedure is provided for research use only. This document is not intended to be used for commercial development or for-profit testing. Use of trade names and commercial sources is for identification only and does not constitute endorsement by the Public Health Service or by the United States Department of Health and Human Services. Please do not distribute this document to other laboratories or commercial entities.**

**1.0 Purpose**

1.1 The purpose of this procedure is to describe clean-up of amplified MRT-PCR product.

**2.0 Definitions**

2.1 EXO I: Exonuclease I

**3.0 Critical Equipment**

3.1 Multichannel Pipettes

3.2 Centrifuge

3.3 Thermocycler

3.4 4°C Refrigerator

3.5 -20°C Freezer

**4.0 Laboratory Materials**

4.1 Aerosol barrier pipette tips

4.2 0.2 ml PCR reaction 8-tube strip and 96-well plate

4.3 Sterile, nuclease free 1.5 ml micro-centrifuge tubes

4.4 Sterile plate sealing foils

4.5 Exonuclease I (20U/μL)

4.6 Nuclease-Free water (not DEPC-treated)

**5.0 Safety Precautions**

5.1 Adhere to the safety guidelines provided in the Biosafety in Microbiological and Biomedical Laboratories and follow all established site-specific safety procedures, including wearing proper personal protective equipment (PPE).

**6.0 Procedure**

6.1 Make Exonuclease (Exo) I master mix by aliquoting the following master mix into a sterile 1.5 micro-centrifuge tube.

Exonuclease I (Exo I) Master Mix	Single reaction	Full 96-well plate
Exonuclease I	0.13 μL	17.2 μL
10X Buffer (supplied in Exo kit)	0.20 μL	27.5 μL
Water	1.67 μL	230.3 μL
<b>Total Volume</b>	<b>2.00 μL</b>	<b>275 μL</b>

6.2 For efficient pipetting, aliquot Exo I master mix evenly into 0.2 mL strip tubes.



- 6.3 Using a multichannel pipette and a new tip for each reaction aliquot 2  $\mu$ L of Exo I master mix into each 25  $\mu$ L MRT-PCR product.
- 6.4 Seal plate securely with foil and briefly spin reaction plate in centrifuge.
- 6.5 Place reaction plate on thermocycler under the “EXO” program with cycling parameters:

Temperature	Time	Repeat
37°C	15min	1x
80°C	15min	1x
10°C	$\infty$	1x

- 6.6 After cycling, the reaction plate can be stored in short-term (overnight to 1 week) at 4°C or long-term (over 1 week to 1 month), store it in -20°C.
- 6.7 Any remaining Exo I master mix can be stored at -20°C.

## 7.0 **Related Procedures**

- 7.1 LP-328 – Multi-segment Reverse Transcription-PCR (MRT-PCR) of Influenza A and B Viruses